

Test Case Planning

Group members	Rin Pereira, Sara Ross, Bobby Khleang, George Williams, Dylan Burton, Omaid Ghafoori
----------------------	--

1. Use Case 1 – *Input medical information*

Test Case ID	TC01 – Valid weight, height, age
Test Objective	Verify that a registered user can enter valid medical information and the system will save it and display a confirmation message.
Preconditions	<ul style="list-style-type: none">• User is a registered user.• User is on the “Medical information” screen (displayed after registration or from User information page).• App/server is available.
Test steps	<ol style="list-style-type: none">1. View the “Medical information” screen showing fields for Weight, Height (ft/in), and Age.2. Enter values for all three fields.3. Click save.
Input Values	<ul style="list-style-type: none">• Weight: 115• Height: 5 ft 5 in• Age: 20
Expected results	<ul style="list-style-type: none">• App verifies inputs against [BR01][BR02].• App stores weight, height, and age.• App displays confirmation text.• On-reopen of the screen (or via profile), the same values are persisted.

Test Case ID	TC02 – Missing field
Test Objective	Verify the app blocks saving when any required field is blank and shows the correct message.
Preconditions	Same as TC01.
Test steps	<ol style="list-style-type: none">1. View the “Medical information” screen.2. Enter only Weight and Age. Leave Height blank.3. Click save.
Input Values	<ul style="list-style-type: none">• Weight: 115• Height: [blank]• Age: 20
Expected results	<ul style="list-style-type: none">• App displays “Height cannot be blank!” (or Weight, Height, Age with correct field name).• No data is saved.• Flow returns to step 2 of the basic path (form remains visible for correction).

Test Case ID	TC03 – Invalid weight
Test Objective	Verify that the app rejects weight values per business rule (must be > 0).

Preconditions	<ul style="list-style-type: none">• Same as TC01.
Test steps	<ol style="list-style-type: none">1. View the “Medication Information” screen.2. Enter 0 for Weight, valid Height and Age.3. Click Save.
Input Values	<ul style="list-style-type: none">• Weight: 0• Height: 5 ft 5 in• Age: 20
Expected results	<ul style="list-style-type: none">• App displays “Invalid weight.”• No data is saved.• Flow returns to step 2 of the basic path.

Test Case ID	TC04 – Invalid age
Test Objective	Verify the app rejects invalid age values per business rules (age > 0 and < 120)
Preconditions	<ul style="list-style-type: none">• Same as TC01.
Test steps	<ol style="list-style-type: none">1. View the “Medical Information” screen.2. Enter valid Weight and Height, set Age to 0.3. Click Save.
Input Values	<ul style="list-style-type: none">• Weight: 115• Height: 5 ft 5 in• Age: 0
Expected results	<ul style="list-style-type: none">• App displays “Invalid Age.”• No data is saved.• Flow returns to step 2 of the basic path.